

Protocol for Dilution Series of Quinine in Solution

Background and Preparation of Stock Solution

For taking accurate measurements of chemicals, a stock solution is prepared with a known concentration. Here, we use mass concentration  $\rho$  defined as

$$\rho = \frac{m}{V}$$

其中m为奎宁的质量，V是溶剂的体积。在下面的例子中，制备具有质量浓度的原溶液

where  $m$  is the mass of quinine and  $V$  is the volume of the solvent. In the following example, a stock solution is prepared with a mass concentration of  $\rho_{stock} = \frac{10 \text{ mg}}{100 \text{ mL}}$  in water.

To quantify the amount of quinine measured in a sample, we make a dilution series.  
为了量化样品中测量的奎宁量，我们做了一个稀释系列

Dissolution Series

To obtain a diluted sample from the stock solution above, we first consider a diluted sample  $\rho_i$  and how we can add additional solvent  $V_i$  to a sample  $V_{stock}$  from the stock solution. This can be accomplished from the following relation

要从上述原液中获得稀释的样品，我们首先考虑稀释的样品，以及如何在样品中添加额外的溶剂从库存溶液

$$\rho_i(V_i + V_{stock}) = \rho_{stock}V_{stock}$$

The overall procedure is to extract a volume  $V_{stock}$  from the stock solution, followed by adding a volume  $V_i$  of the solvent to obtain a desired concentration  $\rho_i$ , we isolate  $V_i$  in the above equation to obtain 整个过程是从原液提取一个 $V_{stock}$ ，然后加入体积 $V_i$ 的溶剂以获得所需浓度，最后分离 $V_i$

$$V_i = \frac{\rho_{stock}V_{stock}}{\rho_i} - V_{stock} = \left(\frac{\rho_{stock}}{\rho_i} - 1\right)V_{stock}$$

We can now list the wanted dilution series of wanted concentrations  $\rho_i$ , compute the volumes we need  $V_i$  when we extract a volume  $V_{stock}$  from the stock solution.

现在可以列出所需浓度的稀释，计算我们需要的体积 $V_i$ ，当我们提取一个体积 $V_{stock}$ 从库存溶液中提取

Table 1: Dilution series of quinine solvated in water at various concentrations

$i$	$\rho_i \left[ \frac{\text{mg}}{100 \text{ mL}} \right]$	$V_{stock} \text{ [mL]}$	$V_i \text{ [mL]}$	$V_{total} = V_{stock} + V_i \text{ [mL]}$
1	7.5	10	3.33	13.33
2	5.0	10	10	20.00
3	2.5	10	30	40.00
4	1.0	10	90	100.0

## 实验方案

### Experimental Protocol

#### 原液

#### Stock Solution

原液是通过测量 $m=10\text{mg}$ 奎宁来配制的。将其转移到试管中，加入 $V=50\text{mL}$ 水使奎宁完全溶解奎宁。此后，加入额外的水，直到试管的最终体积为 $100\text{mL}$ ，得到所需质量浓度的原液。

The stock solution is prepared by measuring  $m = 10 \text{ mg}$  of quinine on a weight scale. Transfer this to a test tube and add  $V = 50 \text{ mL}$  of water to completely dissolve the quinine. Hereafter, additional water is added until the final volume of the test tube reads  $100 \text{ mL}$  yielding a stock solution with the required mass concentration  $\rho_{stock}$ .

#### Dilution Series

For each of the series outlined in Table 1 we perform the following:

- 1) Prepare **test tube**  $i$  for solution with mass concentration  $\rho_i$ .
- 2) Extract  $V_{stock}$  from the stock solution and transfer to test tube  $i$ .
- 3) Extract  $V_i$  from water beaker and transfer to test tube  $i$ .
- 4) Move test tube  $i$  to storage.

1)准备试管，以准备质量浓度为  $\rho_i$  的溶液

2)提取 $V_{stock}$ 从原液中转移到试管中

3)从水杯中提取 $V_i$ ，转移至试管中

4) 将试管移至存储区